

Product Datasheet

Human CKAP4 ELISA Kit (orb404273)

Description	Human Cytoskeleton-associated protein 4(CKAP4) ELISA kit
Reactivity	Human
Range	23.44 pg/mL-1500 pg/mL
Target	CKAP4
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q07065
Sensitivity	5.86 pg/mL
Expiration Date	Please enquire.