

Product Datasheet

Human DNASE1L3 ELISA Kit (orb404289)

Description	Human Deoxyribonuclease gamma(DNASE1L3) ELISA kit
Reactivity	Human
Range	23.44 pg/mL-1500 pg/mL
Target	DNASE1L3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q13609
Sensitivity	5.86 pg/mL
Expiration Date	Please enquire.