

Product Datasheet

Human APEX1 ELISA Kit (orb404323)

Description	Human DNA-(apurinic or apyrimidinic site) lyase (APEX1) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	APEX1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P27695
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.