

Product Datasheet

Human DAG1 ELISA Kit (orb404336)

Description	Human dystroglycan 1 (dystrophin-associated glycoprotein 1) (DAG1) ELISA kit
Reactivity	Human
Range	125 pg/mL-8000 pg/mL
Target	DAG1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q14118
Sensitivity	31.25 pg/mL
Expiration Date	Please enquire.