

Product Datasheet

Human MYLIP ELISA Kit (orb404342)

Description	Human E3 ubiquitin-protein ligase MYLIP(MYLIP) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	MYLIP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q8WY64
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.