

Product Datasheet

Human ENTPD1 ELISA Kit (orb404348)

Description	Human Ectonucleoside triphosphate diphosphohydrolase 1(ENTPD1) ELISA kit
Reactivity	Human
Range	23.44 pg/mL-1500 pg/mL
Target	ENTPD1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P49961
Sensitivity	5.86 pg/mL
Expiration Date	Please enquire.