

Product Datasheet

Human ELTD1 ELISA Kit (orb404353)

Description	Human EGF, latrophilin and seven transmembrane domain-containing protein 1(ELTD1) ELISA kit
Reactivity	Human
Range	25 pg/mL-1600 pg/mL
Target	ELTD1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9HBW9
Sensitivity	6.25 pg/mL
Expiration Date	Please enquire.