

Product Datasheet

Human EDIL3 ELISA Kit (orb404355)

Description	Human EGF-like repeat and discoidin I-like domain-containing protein 3(EDIL3) ELISA kit
Reactivity	Human
Range	125 pg/mL-8000 pg/mL
Target	EDIL3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O43854
Sensitivity	31.25 pg/mL
Expiration Date	Please enquire.