

Product Datasheet

Human ERRFI1 ELISA Kit (orb404396)

Description	Human ERBB receptor feedback inhibitor 1(ERRFI1) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	ERRFI1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9UJM3
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.