

Product Datasheet

Human SULF1 ELISA Kit (orb404407)

Description	Human Extracellular sulfatase Sulf-1(SULF1) ELISA kit
Reactivity	Human
Range	7.8 pg/mL-500pg/mL
Target	SULF1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q8IWU6
Sensitivity	1.95 pg/mL
Expiration Date	Please enquire.