

Product Datasheet

Human FABP5 ELISA Kit (orb404416)

Description	Human Fatty acid-binding protein, epidermal(FABP5) ELISA kit
Reactivity	Human
Range	0.78 ng/mL-50 ng/mL
Target	FABP5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q01469
Sensitivity	0.195 ng/mL
Expiration Date	Please enquire.