

Product Datasheet

Human FGF17 ELISA Kit (orb404426)

Description	Human Fibroblast growth factor 17(FGF17) ELISA kit
Reactivity	Human
Range	3.12 pg/ml-200 pg/ml.
Target	FGF17
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates.
Uniprot ID	O60258
Sensitivity	0.78 pg/ml.
Expiration Date	Please enquire.