

Product Datasheet

Human C4orf7 ELISA Kit (orb404444)

Description	Human Follicular dendritic cell secreted peptide(C4orf7) ELISA kit
Reactivity	Human
Range	7.81 ng/mL-500 ng/mL
Target	C4orf7
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q8NFU4
Sensitivity	1.95 ng/mL
Expiration Date	Please enquire.