

Product Datasheet

Human GABARAP ELISA Kit (orb404464)

Description	Human Gamma-aminobutyric acid receptor-associated protein(GABARAP) ELISA kit
Reactivity	Human
Range	39.06 pg/ml -2500 pg/ml
Target	GABARAP
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	O95166
Sensitivity	9.77 pg/mL
Expiration Date	Please enquire.