

Product Datasheet

Human Glucocorticoid Receptor ELISA Kit (orb404484)

Description	Human Glucocorticoid receptor(NR3C1) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	NR3C1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P04150
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.