

Product Datasheet

Human GRIN1 ELISA Kit (orb404487)

Description	Human Glutamate [NMDA] receptor subunit zeta-1 (GRIN1) ELISA kit
Reactivity	Human
Range	125 pg/mL-8000 pg/mL
Target	GRIN1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q05586
Sensitivity	31.25 pg/mL
Expiration Date	Please enquire.