

## Product Datasheet

### Human QPCT ELISA Kit (orb404491)

|                          |                                                           |
|--------------------------|-----------------------------------------------------------|
| <b>Description</b>       | Human GlutaminyI-peptide cyclotransferase(QPCT) ELISA kit |
| <b>Reactivity</b>        | Human                                                     |
| <b>Range</b>             | 0.312 ng/mL-20 ng/mL                                      |
| <b>Target</b>            | QPCT                                                      |
| <b>Note</b>              | For research use only                                     |
| <b>Application notes</b> | Detection Wavelength: 450 nm                              |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                         |
| <b>Uniprot ID</b>        | <b>Q16769</b>                                             |
| <b>Sensitivity</b>       | 0.078 ng/mL                                               |
| <b>Expiration Date</b>   | Please enquire.                                           |