

Product Datasheet

Human PLA2G12A ELISA Kit (orb404519)

Description	Human Group XIIA secretory phospholipase A2(PLA2G12A) ELISA kit
Reactivity	Human
Range	4.7 pg/mL-300 pg/mL
Target	PLA2G12A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9BZM1
Sensitivity	1.2 pg/mL
Expiration Date	Please enquire.