

## Product Datasheet

### Human GUCA2B ELISA Kit (orb404528)

|                          |                                                                       |
|--------------------------|-----------------------------------------------------------------------|
| <b>Description</b>       | Human guanylate cyclase activator 2B (uroguanylin) (GUCA2B) ELISA kit |
| <b>Reactivity</b>        | Human                                                                 |
| <b>Range</b>             | 31.25 pg/mL-2000 pg/mL                                                |
| <b>Target</b>            | GUCA2B                                                                |
| <b>Note</b>              | For research use only                                                 |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                          |
| <b>Sample Types</b>      | serum, plasma, urine, cell culture supernates                         |
| <b>Uniprot ID</b>        | <b>Q16661</b>                                                         |
| <b>Sensitivity</b>       | 7.8 pg/mL                                                             |
| <b>Expiration Date</b>   | Please enquire.                                                       |