

## Product Datasheet

### Human HSPB2 ELISA Kit (orb404535)

|                          |                                                  |
|--------------------------|--------------------------------------------------|
| <b>Description</b>       | Human Heat shock protein beta-2(HSPB2) ELISA kit |
| <b>Reactivity</b>        | Human                                            |
| <b>Range</b>             | 0.312 ng/mL-20 ng/mL                             |
| <b>Target</b>            | HSPB2                                            |
| <b>Note</b>              | For research use only                            |
| <b>Application notes</b> | Detection Wavelength: 450 nm                     |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                |
| <b>Uniprot ID</b>        | <b>Q16082</b>                                    |
| <b>Sensitivity</b>       | 0.078 ng/mL                                      |
| <b>Expiration Date</b>   | Please enquire.                                  |