

## Product Datasheet

### Human ERVK6 ELISA Kit (orb404553)

|                          |   |
|--------------------------|---|
| <b>Description</b>       | Human HERV-K_7p22.1 provirus ancestral Env polyprotein(ERVK6) ELISA kit |
| <b>Reactivity</b>        | Human   |
| <b>Range</b>             | 0.156 ng/mL-10 ng/mL  |
| <b>Target</b>            | ERVK6   |
| <b>Note</b>              | For research use only   |
| <b>Application notes</b> | Detection Wavelength: 450 nm  |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                                       |
| <b>Uniprot ID</b>        | <a href="#">Q69384</a>  |
| <b>Sensitivity</b>       | 0.039 ng/mL   |
| <b>Expiration Date</b>   | Please enquire.   |