

Product Datasheet

Human ERVK6 ELISA Kit (orb404553)

Description	Human HERV-K_7p22.1 provirus ancestral Env polyprotein(ERVK6) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	ERVK6
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q69384
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.