

Product Datasheet

Human FCGR1A ELISA Kit (orb404557)

Description	Human High affinity immunoglobulin gamma Fc receptor I(FCGR1A) ELISA kit
Reactivity	Human
Range	0.78 ng/mL-50 ng/mL
Target	FCGR1A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	P12314
Sensitivity	0.2 ng/mL
Expiration Date	Please enquire.