

Product Datasheet

Human HAO1 ELISA Kit (orb404590)

Description	Human Hydroxyacid oxidase 1(HAO1) ELISA kit
Reactivity	Human
Range	18.75 pg/mL-1200 pg/mL
Target	HAO1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9UJM8
Sensitivity	4.69 pg/mL
Expiration Date	Please enquire.