

Product Datasheet

Human IGKV4-1 ELISA Kit (orb404594)

Description	Human Ig kappa chain V-IV region(IGKV4-1) ELISA kit
Reactivity	Human
Range	0.156 ng/ml -10 ng/ml
Target	IGKV4-1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates.
Uniprot ID	P06312
Sensitivity	0.039 ng/ml
Expiration Date	Please enquire.