

Product Datasheet

Human PTK7 ELISA Kit (orb404605)

Description	Human Inactive tyrosine-protein kinase 7(PTK7) ELISA kit
Reactivity	Human
Range	12.5 pg/mL-800 pg/mL
Target	PTK7
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q13308
Sensitivity	3.125 pg/mL
Expiration Date	Please enquire.