

Product Datasheet

Human IGFBP6 ELISA Kit (orb404615)

Description	Human Insulin-like growth factor-binding protein 6(IGFBP6) ELISA kit
Reactivity	Human
Range	0.78 ng/mL-50 ng/mL
Target	IGFBP6
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma,tissue homogenates
Uniprot ID	P24592
Sensitivity	0.195 ng/mL
Expiration Date	Please enquire.