

Product Datasheet

Human ITIH5 ELISA Kit (orb404628)

Description	Human Inter-alpha-trypsin inhibitor heavy chain H5(ITIH5) ELISA kit
Reactivity	Human
Range	1.56 ng/mL-100 ng/mL
Target	ITIH5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q86UX2
Sensitivity	0.39 ng/mL
Expiration Date	Please enquire.