

## Product Datasheet

### Human IFI27 ELISA Kit (orb404633)

|                          |                                                                             |
|--------------------------|-----------------------------------------------------------------------------|
| <b>Description</b>       | Human Interferon alpha-inducible protein 27, mitochondrial(IFI27) ELISA kit |
| <b>Reactivity</b>        | Human                                                                       |
| <b>Range</b>             | 28 pg/mL-1800 pg/mL                                                         |
| <b>Target</b>            | IFI27                                                                       |
| <b>Note</b>              | For research use only                                                       |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                                |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates                                           |
| <b>Uniprot ID</b>        | <b>P40305</b>                                                               |
| <b>Sensitivity</b>       | 7 pg/mL                                                                     |
| <b>Expiration Date</b>   | Please enquire.                                                             |