

## Product Datasheet

### Human IFNE ELISA Kit (orb404635)

|                          |                                             |
|--------------------------|---------------------------------------------|
| <b>Description</b>       | Human Interferon epsilon(IFNE)<br>ELISA kit |
| <b>Reactivity</b>        | Human                                       |
| <b>Range</b>             | 15.6 pg/mL-1000 pg/mL                       |
| <b>Target</b>            | IFNE                                        |
| <b>Note</b>              | For research use only                       |
| <b>Application notes</b> | Detection Wavelength: 450 nm                |
| <b>Sample Types</b>      | serum, plasma, tissue<br>homogenates        |
| <b>Uniprot ID</b>        | <b>Q86WN2</b>                               |
| <b>Sensitivity</b>       | 3.9 pg/mL                                   |
| <b>Expiration Date</b>   | Please enquire.                             |