

Product Datasheet

Human IFIT1 ELISA Kit (orb404645)

Description	Human Interferon-induced protein with tetratricopeptide repeats 1(IFIT1) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	IFIT1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P09914
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.