

Product Datasheet

Human ILF3 ELISA Kit (orb404646)

Description	Human Interleukin enhancer-binding factor 3(ILF3) ELISA kit
Reactivity	Human
Range	15.6 pg/mL-1000 pg/mL
Target	ILF3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma,tissue homogenates, cell lysates
Uniprot ID	Q12906
Sensitivity	3.9 pg/mL
Expiration Date	Please enquire.