

Product Datasheet

Human IL22RA2 ELISA Kit (orb404664)

Description	Human Interleukin-22 receptor subunit alpha-2(IL22RA2) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	IL22RA2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q969J5
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.