

## Product Datasheet

### Human IL27RA ELISA Kit (orb404665)

|                          |                                                               |
|--------------------------|---------------------------------------------------------------|
| <b>Description</b>       | Human Interleukin-27 receptor subunit alpha(IL27RA) ELISA kit |
| <b>Reactivity</b>        | Human                                                         |
| <b>Range</b>             | Request Information                                           |
| <b>Target</b>            | IL27RA                                                        |
| <b>Note</b>              | For research use only                                         |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                  |
| <b>Uniprot ID</b>        | <b>Q6UWB1</b>                                                 |
| <b>Sensitivity</b>       | Request Information                                           |
| <b>Expiration Date</b>   | Please enquire.                                               |