

Product Datasheet

Human IL1F5 ELISA Kit (orb404673)

Description	Human Interleukin-36 receptor antagonist protein(IL1F5) ELISA kit
Reactivity	Human
Range	31.25 pg/ml - 2000 pg/ml
Target	IL1F5
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q9UBH0
Sensitivity	7.8 pg/ml
Expiration Date	Please enquire.