

Product Datasheet

Human KEAP1 ELISA Kit (orb404686)

Description	Human Kelch-like ECH-associated protein 1(KEAP1) ELISA kit
Reactivity	Human
Range	25 pg/mL-1600 pg/mL
Target	KEAP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q14145
Sensitivity	6.25 pg/mL
Expiration Date	Please enquire.