

Product Datasheet

Human KRT6A ELISA Kit (orb404689)

Description	Human Keratin, type II cytoskeletal 6A(KRT6A) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	KRT6A
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P02538
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.