

Product Datasheet

Human LTBP1 ELISA Kit (orb404708)

Description	Human Latent-transforming growth factor beta-binding protein 1(LTBP1) ELISA kit
Reactivity	Human
Range	0.156 ng/mL-10 ng/mL
Target	LTBP1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	Q14766
Sensitivity	0.039 ng/mL
Expiration Date	Please enquire.