

Product Datasheet

Human LRRC17 ELISA Kit (orb404725)

Description	Human Leucine-rich repeat-containing protein 17(LRRC17) ELISA kit
Reactivity	Human
Range	23.44 pg/mL-1500 pg/mL
Target	LRRC17
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q8N6Y2
Sensitivity	5.86 pg/mL
Expiration Date	Please enquire.