

Product Datasheet

Human LECT2 ELISA Kit (orb404731)

Description	Human Leukocyte cell-derived chemotaxin-2(LECT2) ELISA kit
Reactivity	Human
Range	1.56 ng/mL-100 ng/mL
Target	LECT2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell culture supernates
Uniprot ID	O14960
Sensitivity	0.39 ng/mL
Expiration Date	Please enquire.