

Product Datasheet

Human SERPINB1 ELISA Kit (orb404732)

Description	Human Leukocyte elastase inhibitor(SERPINB1) ELISA kit
Reactivity	Human
Range	0.625 ng/mL-40 ng/mL
Target	SERPINB1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P30740
Sensitivity	0.156 ng/mL
Expiration Date	Please enquire.