

Product Datasheet

Human FCGR3A ELISA Kit (orb404761)

Description Human Low affinity immunoglobulin gamma Fc region receptor III-A(FCGR3A) ELISA kit

Reactivity Human

Range 47 pg/mL-3000 pg/mL

Target FCGR3A

Note For research use only

Application notes Detection Wavelength: 450 nm

Sample Types serum, plasma, tissue homogenates

Uniprot ID **P08637**

Sensitivity 11.75 pg/mL

Expiration Date Please enquire.