

Product Datasheet

Human LYPD1 ELISA Kit (orb404767)

Description	Human Ly6/PLAUR domain-containing protein 1(LYPD1) ELISA kit
Reactivity	Human
Range	25 pg/mL-1600 pg/mL
Target	LYPD1
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q8N2G4
Sensitivity	6.25 pg/mL
Expiration Date	Please enquire.