

Product Datasheet

Human LYPD3 ELISA Kit (orb404768)

Description	Human Ly6/PLAUR domain-containing protein 3(LYPD3) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	LYPD3
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates
Uniprot ID	O95274
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.