

Product Datasheet

Human LIPA ELISA Kit (orb404777)

Description	Human Lysosomal acid lipase/cholesteryl ester hydrolase(LIPA) ELISA kit
Reactivity	Human
Range	62.5 pg/mL-4000 pg/mL
Target	LIPA
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P38571
Sensitivity	15.6 pg/mL
Expiration Date	Please enquire.