

Product Datasheet

Human LOXL2 ELISA Kit (orb404784)

Description	Human Lysyl oxidase homolog 2(LOXL2) ELISA kit
Reactivity	Human
Range	37.5 pg/mL-2400 pg/mL
Target	LOXL2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9Y4K0
Sensitivity	9.375 pg/mL
Expiration Date	Please enquire.