

## Product Datasheet

### Human MATR3 ELISA Kit (orb404799)

|                          |                                                   |
|--------------------------|---------------------------------------------------|
| <b>Description</b>       | Human Matrin-3(MATR3) ELISA kit                   |
| <b>Reactivity</b>        | Human                                             |
| <b>Range</b>             | 23.5 pg/mL-1500 pg/mL                             |
| <b>Target</b>            | MATR3                                             |
| <b>Note</b>              | For research use only                             |
| <b>Application notes</b> | Detection Wavelength: 450 nm                      |
| <b>Sample Types</b>      | serum, plasma,tissue<br>homogenates, cell lysates |
| <b>Uniprot ID</b>        | <b>P43243</b>                                     |
| <b>Sensitivity</b>       | 5.8 pg/mL                                         |
| <b>Expiration Date</b>   | Please enquire.                                   |