

Product Datasheet

Human MAGED2 ELISA Kit (orb404817)

Description	Human Melanoma-associated antigen D2(MAGED2) ELISA kit
Reactivity	Human
Range	31.25 pg/mL-2000 pg/mL
Target	MAGED2
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	Q9UNF1
Sensitivity	7.81 pg/mL
Expiration Date	Please enquire.