

## Product Datasheet

### Human DNAJC19 ELISA Kit (orb404841)

|                          |                                                                                        |
|--------------------------|----------------------------------------------------------------------------------------|
| <b>Description</b>       | Human Mitochondrial import inner membrane translocase subunit TIM14(DNAJC19) ELISA kit |
| <b>Reactivity</b>        | Human                                                                                  |
| <b>Range</b>             | 28 pg/mL-1800 pg/mL                                                                    |
| <b>Target</b>            | DNAJC19                                                                                |
| <b>Note</b>              | For research use only                                                                  |
| <b>Application notes</b> | Detection Wavelength: 450 nm                                                           |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates, cell lysates                                        |
| <b>Uniprot ID</b>        | <b>Q96DA6</b>                                                                          |
| <b>Sensitivity</b>       | 7 pg/mL                                                                                |
| <b>Expiration Date</b>   | Please enquire.                                                                        |