

Product Datasheet

Human MT-ND6 ELISA Kit (orb404895)

Description	Human NADH-ubiquinone oxidoreductase chain 6(MT-ND6) ELISA kit
Reactivity	Human
Range	23.44 pg/mL-1500 pg/mL
Target	MT-ND6
Note	For research use only
Application notes	Detection Wavelength: 450 nm
Sample Types	serum, plasma, tissue homogenates, cell lysates
Uniprot ID	P03923
Sensitivity	5.86 pg/mL
Expiration Date	Please enquire.